


Fisher & Paykel Healthcare

Mask Family Seal Size Measurements

Nasal Masks



F&P Eson 2




F&P Eson

Nose Width

Less than 3.7 cm (1.46 inches)	Between 3.7 cm and 4.1 cm (1.46 inches and 1.61 inches)	Greater than 4.1 cm (1.61 inches)
-----------------------------------	--	--------------------------------------

Nose Height	Less than 4.4 cm (1.73 inches)	S	M	L
	Between 4.4 cm and 5.2 cm (1.73 inches and 2.05 inches)	M	M	L
	Greater than 5.2 cm (2.95 inches)	L	L	L

Full Face Masks



F&P Vitera



F&P Simplus

Small	Less than 8.9 cm (3.5 inches)
Medium	Between 8.9 cm and 10 cm (3.5 inches and 3.94 inches)
Large	Greater than 10 cm (3.94 inches)

Measuring these points for sizing in a mirror may make the process easier for you.

These products may not be right for you. Read the warnings and instructions for use before purchase.

For further instructions on sizing, please refer to the Eson™ 2 and Simplus™ sizing guide REF 185048403 or Vitera™ sizing guide REF 612027

F&P, Eson, Vitera and Simplus are trademarks of Fisher & Paykel Healthcare Limited.
For patent information, see www.fphcare.com/ip

Prescription only. Use only as directed. Always follow the instructions for use.
Your healthcare professional will advise you whether this product is suitable for you/your condition.